

Installation Guide

CTS FM50 Stack Mount

The purpose of this guide is to:


- Describe attaching Stack Mount to CTS (page 1-2).
- Describe attaching CTS Stack Mount to FM50 (page 2).
- Describe Routine Maintenance (page 3).



Parts Reference

The following parts are included in this installation kit :

| Item # | Description | Qty |
|--------|--|-----|
| 1 | Adapter Plate | 1 |
| 2 | Base Plate | 1 |
| 3 | M4 x 10mm Flat Head Machine Screw (FHMS) | 4 |
| 4 | M3.5 X 12MM (PHMS) | 6 |




Tools Required:


- Phillips Screwdriver

Attaching Stack Mount to CTS

1. Using a Phillips screwdriver, attach Adapter Plate to the bottom of the CTS with (4) M3.5 X 12mm PHMS.




2. Attach Base Plate to Adapter plate with (4) M4 x 10mm FHMS.



Attaching CTS Stack Mount to FM50

Installation Note: Plastic inserts may need to be removed from the (2) threaded mounting holes on top of the FM50.

1. Secure CTS Stack Mount to FM50 with (2) M3.5 x 12mm PHMS as shown below.



Routine Maintenance

Periodically inspect all fasteners associated with the mounting system. Tighten as necessary for optimal operation and safety.

Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.
3. Do not use steel wool or other abrasive material to clean the mounting assembly.
4. Damage caused by the use of unapproved substances or processes will not be covered by warranty. We recommended that you test any cleaning solution on a small area of the mounting assembly that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the mounting assemblies. Wipe any cleaning agents off the mounting assembly immediately using a water-dampened cloth. Dry all mounting assemblies thoroughly after cleaning.

CAUTION: GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted instruments or accessory equipment, refer to the specific instructions delivered with those products.