

### Installation Guide

Datex-Ohmeda Aespire Core Arm Mounting Assembly for Bottom Mounted Devices  
 For use with DX-0033-01 Ventilator Display Arm - Refer to Qualified Personnel




25lbs/11.3kg Maximum Device Weight

### Parts Reference



The following parts and hardware are included with this installation kit (see photos for parts, hardware not shown):

Item #	Description	Qty
1	Mounting Plate/Riser Assembly	1
2	#10-32 x 2" Socket Head Cap Screw (SHCS)	3
3	5/32" Hex Key	1



**Tools Required:** 5/32" hex key (provided)

1. Identify three (3) mounting holes on top of the Ventilator Display Arm (left side of Aespire) shown below. These mounting holes will be located directly behind the ventilator display.
2. Align the threaded mounting holes in the end of the Riser with the three (3) mounting holes in the Ventilator Display Arm. Using the 5/32" [4-mm] hex key, secure the Mounting Plate/Riser Assembly with three (3) #10-32 x 2" SHCS inserted from the bottom of the Arm. **Caution:** Do not over-tighten these screws.




## Mounting Device on Slide-In Type Mounting Plate

1. Pull Spring Plunger at the front of the Mounting Plate and slide device laterally into the Mounting Plate until the Plunger snaps into the clearance hole at the front of the Mounting Adapter (**Note:** Mounting Adapter, if needed, sold separately).
2. Tighten the Nylon Locking Screws located on bottom or rear of the Mounting Plate

## Adjusting Tilt Tension on the Mounting Plate/Riser Assembly

Adjust tilt after device is installed. Two (2) Tilt Adjustment screws are located on the front of the Mounting Assembly. Using a 5/32" [4-mm] hex key (provided), equally loosen or tighten the socket head screws until the monitor can be tilted manually *and* maintains its position after tilting.



## Routine Maintenance

Periodically inspect all fasteners associated with the mounting system. Tighten or adjust fasteners as necessary for optimal operation and safety.

## Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone or trichloroethylene.
3. Do not use Steel wool or other abrasive material.
4. Damage caused by the use of unapproved substances or processes will not be warranted. It is recommended that you test any cleaning solution on a small area of the mounting assembly that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the mounting assembly. Wipe any cleaning agents off of the assembly immediately, using water dampened cloth. Dry the assembly thoroughly after cleaning.

**CAUTION:** GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted devices or accessory equipment, refer to the specific instructions delivered with those products.