

Installation Guide

Philips FMS and X2 Flush Channel Mount (Refer to qualified personnel)

The purpose of this guide is to:

1. Describe attachment of Slide Adapter to FMS/X2 unit (page 2).
2. Describe mounting of flush mount to the wall channel (page 2).




Parts Reference

The following parts and hardware are included in this installation kit:

Item #	Description	Qty
1	Slide Extension Adapter	1
2	M4 x 12mm Pan Head Machined Screw (PHMS)	2
3	1/8" Hex Wrench	1



1



Tools required

- 1/8" hex wrench (provided)
- Phillips Screwdriver (not provided)



WARNING: Do not position support arm or device above a patient.



Before mounting the arm:

1. Verify that the Channel has been installed and approved in accordance with the GCX Wall Channel Installation Instructions.
2. This slide adapter is rated for 10 lbs/4.5 kg. It is not recommended that you attempt to use the arm for weights outside of this range.
3. If you are unsure of your application, please contact a GCX product specialist at (800) 228-2555 for assistance.

Failure to periodically inspect and adjust these functions may result in damage to equipment or injury.


Attaching Slide Extension Adapter to FMS/X2

Partially install M4 x 12mm PHMS (2) by hand, and give space to allow Slide Extension Adapter to slide in. Insert and slide adapter as pictured, and tighten the screws with a Philips head screwdriver until adapter is secure to the FMS/X2.



Mounting FMS/X2 to Channel

Slide the adapter from the top of the channel down to the desired position and tighten screws with the hex wrench until stable.



Routine Maintenance

Periodically inspect all screws associated with the mounting system. Tighten or adjust screws as necessary for optimal operation and safety.

Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone or trichloroethylene.
3. Do not use steel wool or other abrasive material.
4. Damage caused by the use of unapproved substances or processes will not be warranted. It is recommended that you test any cleaning solution on a small area of the mounting assembly that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the mounting assembly. Wipe any cleaning agents off of the mounting assembly immediately, using water dampened cloth. Dry the mounting assembly thoroughly after cleaning.

CAUTION: GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted devices or accessory equipment, refer to the specific instructions delivered with those products before installing the arm in the Channel.