

Installation Guide
Tympanic Thermometer 1.25" Diameter Post Mount



Parts Reference

The following parts and hardware are included with this installation kit (hardware not shown):

Item #	Description	Qty
1	Front Clamping Block	1
2	Rear Clamping Block	1
3	M4 x 12mm Pan Head Machine Screw (PHMS)	2
4	#10-32 x 9/16" PHMS	2



Front Clamping Block




Back Clamping Block

Tools Required: Phillips screwdriver (not provided).


Mounting the Tympanic Thermometer Base Station

Installation Note: Install and route Base Station cables prior to fastening Clamping Block to Base Station.

- Using a Phillips screwdriver, fasten the front Clamping Block to Base Station with two (2) M4 x 12 mm PHMS as shown below.



- Using a Phillips screwdriver, secure the Base Station mounting assembly to post with the back Clamping Block and two (2) 10-32 x 9/16" PHMS as shown below.



Routine Maintenance

Periodically inspect all fasteners associated with the mounting system. Tighten or adjust fasteners as necessary for optimal operation and safety.

Cleaning the Mounting Assembly

- The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
- The surface finish will be permanently damaged by strong chemicals and solvents such as acetone or trichloroethylene.
- Do not use Steel wool or other abrasive material.
- Damage caused by the use of unapproved substances or processes will not be warranted. It is recommended that you test any cleaning solution on a small area of the arm that is not visible to verify compatibility.
- Never submerge or allow liquids to enter the mounting assembly. Wipe any cleaning agents off of the arm immediately, using water dampened cloth. Dry the assembly thoroughly after cleaning.

CAUTION: GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted devices or accessory equipment, refer to the specific instructions delivered with those products before installing the arm in the Channel.