

Installation Guide

Spacelabs 91393 Adapter Plate

The purpose of this guide is to:

- Describe attachment of Adapter Plate to the Spacelabs 91393 (page 1).
- Describe mounting Spacelabs 91393 onto Slide-on Mounting Plate (page 2).
- Describe Routine Maintenance (page 3).

Parts Reference

The following parts are included in this installation kit (hardware not shown):

| Item # | Description | Qty |
|--------|---|-----|
| 1 | Adapter Plate | 1 |
| 2 | 8-32 x 1.25" Flat Head Machine Screw (FHMS) | 4 |




Tools Required:

- Phillips screwdriver (not provided)


Attaching Adapter Plate to Spacelabs 91393

Using a Phillips screwdriver, attach Adapter Plate to the bottom of the Spacelabs 91393 with (4) 8-32 x 1.25" FHMS as shown.



Mounting Spacelabs 91393 on Slide-on Mounting Plate

Slide Mounting Adapter into Mounting Plate while pulling spring-loaded plunger at bottom of Plate. Plunger will snap into clearance hole when the 91393 is correctly mounted.



Routine Maintenance

Periodically inspect all fasteners associated with the mounting system. Tighten as necessary for optimal operation and safety.

Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.
3. Do not use steel wool or other abrasive material to clean the mounting assembly.
4. Damage caused by the use of unapproved substances or processes will not be covered by warranty. We recommended that you test any cleaning solution on a small area of the mounting assembly that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the mounting assemblies. Wipe any cleaning agents off the mounting assembly immediately using a water-dampened cloth. Dry all mounting assemblies thoroughly after cleaning.

CAUTION: GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted instruments or accessory equipment, refer to the specific instructions delivered with those products.