

VHM-Extension Retrofit Kit



- VHM can be mounted above or below the extension
- Folds flush against the wall for minimal foot print while stowed
- Dual pivoting extension is ideal for obtaining additional reach and allows fold over feature
- Optional rotation stop at front pivot of extension +/-100 degrees (200 degrees total range)
- Cable management

8"/20.3 cm VHM-Extension Kit

WS-0005-74

14"/35.6 cm VHM-Extension Kit

WS-0005-75



Fold-over feature

3875 Cypress Drive
 Petaluma, California 94954-5635
 United States
 800-228-2555
 707-773-1100 outside US
 707-773-1180 fax

Postbus 10462
 5000 JL Tilburg
 Netherlands
 +31 (0) 88 627 2600
 +31 (0) 88 627 2609 fax

No. 7, Lane 116, Wugong 3rd Road
 Wugu Township, Taipei County 248
 Taiwan (R.O.C.)
 (886) 2 2298 2842
 (886) 2 8990 3052 fax

gcx.com
 one tool. one source.