


Light Duty Roll Stand Package for Bottom Mounted Devices



GCX P/N: RS-0025-02C

- ❑ Load capacity up to 12 lbs (5.5 kg)*
- ❑ Overall Height 44" (111.8 cm)
- ❑ 21" (53.3 cm) diameter aluminum base
- ❑ 3" (7.6 cm) single wheel casters (2 locking, 3 non-locking) for easy transition over uneven surfaces
- ❑ 5" (12.7 cm) slide-in type mounting plate
- ❑ Mounting plate fixed at 15 degree angle
- ❑ 6" (15.2 cm) basket (12" w x 8" d x 6" h) with hooks for cable management
- ❑ Handle with padded grip

*Counterweight options available for heavier monitors

Contact GCX for mounting adapter solutions.

GCX® offers a complete range of mounting solutions. Please contact us or visit www.gcx.com for additional information.

GCX
INSTRUMENT MOUNTING SYSTEMS

3875 Cypress Drive
Petaluma, CA 94954-5635
707.773.1100
800.228.2555 US/CAN
Fax: 707.773.1180
E-mail: sales@gcx.com
www.gcx.com

GCX is an ISO9001 registered company.



GCX, PolyMount and PolyQuip are registered trademarks of GCX Corporation.

© Copyright 2005

2/05-1



Go to www.gcx.com for the latest information about our products and services.