



GCX Fetal Monitoring Workstations

A full range of solutions from the experts in wall mounts and carts for clinical and IT devices

Up to
\$131,754

excess cost per HAI incident¹

Medical-grade cleanability for superior infection control. Built with metal, wood and minimal or no thermofoil finishes—GCX Fetal Monitoring Workstations are designed to meet Joint Commission guidance for equipment and disinfectant selection.

¹ *Open Forum Infectious Diseases, Volume 9, Issue Supplement 2, December 2022*

Broad portfolio to meet every space and budget requirement

Choose right-sized, right-priced, fully-configurable solutions for Labor & Delivery suites, OR, ER, triage and other hospital locations and applications.

Built to last from medical-grade materials

All GCX Fetal Monitoring Workstations are engineered for years of dependable performance in the most challenging hospital environments, and backed by industry-standard warranties.

Supports GE, Philips and other leading manufacturers

All workstations are compatible with fetal monitors from GE, OBIX, Philips and other leading brands, including many monitors purpose-built for other carts.



Jaco 520 Series Compact Fetal Monitoring Carts *Space-saving, cost-effective fetal monitor mobility*

- Our smallest footprint (16" x 16") and lowest price point
- Built with aircraft-grade metal; high-gloss antimicrobial finish
- Two (2) compact storage drawers and large worksurface
- Height adjustable work surface (41" to 48")
- Fits easily in hall and out of the way when not in use
- Choice of single or dual monitor stands with tilt and pan LDC display adjustments



GCX RC Series Classic Fetal Monitoring Carts

Trusted mobile fetal monitoring solution for OR, ER and triage

- Compact 18" x 18" footprint
- Up to 2 VESA display mounts with tilt, pan and landscape-to-portrait adjustments
- 2-3 spacious storage drawers, printer paper catch and 4 quick-release baskets
- Integrated channel accepts most GCX channel mounting components
- Variable height keyboard mount moves easily out of the way
- Locking 5" casters for smooth movement over all flooring types and transitions



Classic Fetal Monitoring Wall Mount

Full featured monitoring workstation—ideal when floor space is limited

- Mounts on trusted GCX wall mounting channel; accepts all GCX mounting components
- Several configurations available to meet a variety of height adjustment requirements
- Up to 2 VESA display mounts with tilt, pan and landscape-to-portrait adjustments
- 2 spacious storage drawers, printer paper catch and 4 quick-release baskets
- Variable height keyboard mount moves easily out of the way

GCX MC Series

Mobile Fetal Monitoring Workstation

Comprehensive mobile fetal monitoring workstation for Labor and Delivery (L&D) suites

- Designed and refined in consultation with L&D nurses
- Up to 2 VESA display mounts with tilt, pan and landscape-to-portrait adjustments
- Ability to mount 2 displays, side-by-side or stacked
- Large worksurface with handle, integrated gel bottle holders, and spill channel
- Up to 6 spacious storage drawers with printer paper catch
- Choice of pull-out or sit-to-stand height adjustable side mounted keyboard
- Comforting patient-friendly touches: soft nighttime charting light, rounded corners, and available wood-tone accents
- Locking 5" casters for smooth movement over all flooring types and transitions



Product comparison

| | GCX Jaco UL 520 FM Cart | GCX RC Series FM Cart | GCX FM Wall Mount | GCX MC Series FM Cart |
|---------------------------------|--|---|---|--|
| LCD (VESA) Mounts | 1 or 2 (Side-by-Side) | 1 or 2 (Side-by-Side or Stacked) | 1 or 2 (Side-by-Side or Stacked) | 1 or 2 (side-by-side, or Stacked with Optional Height Adjustment) |
| CPU Mount | Yes | Yes | Yes | Yes |
| Max. Drawers | 2 | 2 | 2 | 6* |
| Paper Catch | Yes | Yes | Yes | Yes |
| Keyboard/Mouse Tray | Pull-Out | Height-Adjustable | Optional Height-Adjustable | Pull-Out or Side Mounted/Height-Adjustable |
| Worksurface Height Adjustment | 7" | N/A | N/A | N/A |
| Keyboard Tray Height Adjustment | 7" | 12" | 12" | 12" |
| Lighting | Keyboard Light (Powered Carts Only) | N/A | N/A | Keyboard/Worksurface Light |
| Cable Management | Yes | Yes | Yes | Yes |
| Lead Management | N/A | Optional | Optional | Yes |
| Footprint | 16"W x 16"D | 18"W x 18"D | N/A | 25"W x 25"D |
| Color | White | White | White | White |
| Accent Color Options | 10 Standard Handle Colors | N/A | N/A | Gray or 3 Woodtones (Worksurface) |
| Accessory Mounts | Barcode Scanners, Baskets, Wipes Holder, Utility Cup, Printer Mounts | Multiple Telemetry, Barcode Scanner and Basket Mounts | Multiple Telemetry, Barcode Scanner and Basket Mounts | Multiple Telemetry, Barcode Scanner and Basket Mounts; Gel Bottle Holdster |

*With Side-Mounted Keyboard Tray. Max. 3 Drawers with Pull-Out Keyboard Tray.



sales@gcx.com

3875 Cypress Drive
Petaluma, California 94954
800.228.2555

140 Constitution Blvd.
Franklin, MA 02038
800.649.2278