

## Philips PM6000 Series on M Series Pivot Arm Channel Mount

### FEATURES

- Offers a pivot point at the wall channel for lateral adjustability
- Features range of motion from flush to full extension, and cable management
- Tilt mechanism
- Accessory bracket for Evoline rack

### PRODUCT DETAIL

- 8" / 20.3 cm M Series Pivot Arm with Slide-In Mounting Plate
- Adapter plate for device included
- Tilt mechanism maximum load: 30 lb / 13.6 kg
- Total arm maximum load: 60" / 27.2 kg



GCX is an ISO 13485 registered company (TUV-USA). GCX, PolyMount, and PolyQuip are registered trademarks, and VHRC, VHC, VHM, VHRS, PRC, M Series, and Stor-Lock are trademarks, of GCX Corporation. All other brand names, product names, or trademarks belong to their respective holders. 2025 GCX Corporation

[www.gcx.com](http://www.gcx.com)

Find this product at <https://www.gcx.com/product/philips-evoline-pm6000-on-m-series-arm/>