



# Fleet Management Software for EVO Medical Computer Carts

Put your fleet data to work with real-time cart insights and robust integration capabilities

METRIQ by GCX FLEET MANAGEMENT DASHBOARD										
Keyword Search		Access Level	Status		Group		CLEAR			
Unique ID	Status	Feature	Cart Name	Description	Group	Battery 1	Battery 2	Last Reported Location	Data Timestamp	Actions
GSD-CC22B811AA07	Idle	L500	FSC-ER-F2-E-04	ER-F2	ER-F2	15 %	52 %	1st Floor, Emergency Room	02/25/2025 07:30:25	
GSD-CC22B811AA05	Idle	L500	FSC-ER-F2-E-02	ER-F2	ER-F2	65 %	11 %	1st Floor, Emergency Room	02/25/2025 07:30:25	
GSD-CC22B811AA04	Idle	HS	FSC-ER-F2-E-01	ER-F2	ER-F2	80 %	44 %	1st Floor, Emergency Room	02/25/2025 07:30:25	
GSD-CC22B811AA01	Idle	HS	FSC-ER-F1-W-02	ER-F1	ER-F1	26 %	2 %	1st Floor, Emergency Room	02/25/2025 07:30:25	
GSD-CC22B811AA03	Idle	L500	FSC-ER-F1-W-04	ER-F1	ER-F1	66 %	94 %	1st Floor, Emergency Room	02/25/2025 07:30:25	
GSD-CC22B811AA00	Idle	HS	FSC-ER-F1-W-01	ER-F1	ER-F1	7 %	42 %	1st Floor, Emergency Room	02/25/2025 07:30:25	
GSD-CC22B811AA02	Active	HS	FSC-ER-F1-W-03	ER-F1	ER-F1	17 %	65 %	1st Floor, Emergency Room	02/25/2025 07:30:25	
GSD-CC22B811AA09	Idle	L250	FSC-LND-F1-N-03	LND-F1	LND-F1	24 %	69 %	3rd floor, L&D	02/25/2025 07:30:25	
GSD-CC22B811AA08	Active	L250	FSC-LND-F1-N-02	LND-F1	LND-F1	75 %	81 %	3rd floor, L&D	02/25/2025 07:30:25	
GSD-CC22B811AA13	Idle	L250	FSC-LND-F1-N-06	LND-F1	LND-F1	1 %	24 %	3rd floor, L&D	02/25/2025 07:30:25	

## Real-time cart data and analytics

METRIQ's real-time cart location, access and utilization data help you track down drifted carts, understand usage patterns, redeploy carts for optimal utilization and right-size your fleet. A predictive battery life algorithm, easy access to warranty data and other tracking features help you plan and budget accurately.

## Microsoft Active Directory integration

METRIQ uses your existing Active Directory structure to manage access levels, assign roles and distinguish application users from cart users. Authenticate users with your existing single sign-on (and without additional passwords to remember and manage).

## Seamless service ticketing integration

Use out-of-the-box with your ServiceNow ticketing and with GCX tech support. Cart users can initiate detailed service requests with the touch of a button; IT staff can quickly escalate service requests to GCX when needed. Automate maintenance scheduling and reminders to maximize cart life and performance.

## Fleet data API

Extend the value of your cart fleet data by integrating METRIQ with asset management, preventative maintenance/CMMS, healthcare facilities management (HFM) and other mission critical systems.

## Cloud-based convenience and security

With METRIQ there's no software to install on client devices—just log in from any browser to access your fleet data. Keep fleet management and EHR separate for streamlined HIPAA compliance. And protect your data with Amazon Web Services KMS, API Gateway, GuardDuty and CloudWatch.

[sales@gcx.com](mailto:sales@gcx.com)

3875 Cypress Drive  
Petaluma, California 94954  
800.228.2555

140 Constitution Blvd.  
Franklin, MA 02038  
800.649.2278